

## Tyvek® Boot Cover

### Extended foot protection

Cover extends over your ankle for extra protection. Elastic opening makes cover easy to put on and take off.



86225-52

Catalog number	Description	Size	Qty/cs	Price/cs
<a href="#">K-86225-52</a>	Tyvek boot covers	Universal	200	

## Hazmat Boot Covers

### Ideal for hazardous materials handling

Hazmat boot covers are made from 0.70-mm latex rubber and are 100% waterproof—ideal for hazardous material handling applications. Ribbed and textured soles help prevent slipping. Covers are 12"H.



86235-51

Catalog number	Description	Size	Price per pair	
			1 pr	12 pr
<a href="#">K-86235-51</a> <a href="#">K-86235-52</a>	Hazmat boot covers	Large (7-9) X-large (10-15)		

## Acrylic Shoe Cover Dispenser

### Keep your work area tidy

Clear acrylic dispenser holds over 50 pairs of disposable shoe covers. Cutout opening allows convenient access at front—easy fill flap at top. Dispenser can be set on the benchtop or wall mounted. Measures 12"W x 4"H x 8"D.



86456-10

Catalog number	Capacity	Dimensions	Price
<a href="#">K-86456-10</a>	50 pairs	12"W x 14"H x 8"D	

## Batta Rubber Overboots

### Wear over shoes for added protection

Shoes and boots feature a four-way cleated sole that gives you superior slip resistance. Reinforced toe and heel help boots last longer. Ribs on back of boot quickly identify size. Made in USA.



86235-21

86235-31

Catalog number	Description	Size	Height	Price per pair	
				1 pr	6 pr
<a href="#">K-86235-20</a> <a href="#">K-86235-21</a> <a href="#">K-86235-22</a>	Overshoe	M (8-9) L (10-11) XL (12-13)	4"		
<a href="#">K-86235-30</a> <a href="#">K-86235-31</a> <a href="#">K-86235-32</a>	Overboot	M (8-9) L (10-11) XL (12-13)	10"		

## Polypropylene (PP) Shoe Covers

### Contain particulates from shoe soles

These shoe covers are made of 100% spunbond PP. The non-woven fabric is available in two styles: regular or skid-free sole. Both styles come in extra large and have an elastic band at the ankles to hold them to your feet. Shoe covers with skid-free soles provide greater traction for improved safety.

Regular soles  
86227-82



Skid-free soles  
86227-92



Catalog number	Description	Size	Qty/cs	Price per case	
				1 cs	12 cs
<a href="#">K-86227-82</a> <a href="#">K-86227-92</a>	Regular sole Skid-free sole	Extra large	300		

## Polyethylene Shoe Covers

### Economical splash protection

These shoe covers are an economical alternative when a low particulate material is required to protect against splashes. Covers are made of lint-free, water-resistant, low-density polyethylene.



86228-20

Catalog number	Description	Size	Qty/cs	Price per case	
				1 cs	12 cs
<a href="#">K-86228-20</a>	Polyethylene shoe covers	Universal	1000		

## Polylatex Shoe Covers

### Liquid resistant coverage

These shoe covers are made of durable, water resistant, lint-free blue polylatex. Polylatex offers better protection from liquid and more durability than other materials. Available in universal and extra-large sizes.



86228-10

Catalog number	Description	Size	Qty/cs	Price per case	
				1 cs	12 cs
<a href="#">K-86228-10</a> <a href="#">K-86228-12</a>	Polylatex shoe covers	Universal Extra large	300		

## Tyvek® Shoe Cover

### Rugged protection from punctures and tears

Use to keep shoes covered and protected. Elastic opening at the ankle allows quick and easy removal.



86225-50

Catalog number	Description	Size	Qty/cs	Price/cs
<a href="#">K-86225-50</a>	Tyvek shoe covers	Universal	100	